

SafeTEE

Our *PROVEN* 4-step solution for mitigating costly TEE probe failures



Supporting TEE probes can take a bite out of any budget. But, did you know that most catastrophic TEE probe failures are fully avoidable?

Our clinical and technical experts provide you with the tools and expertise to bite back.

SafeTEE is a 4-step program proven to significantly reduce TEE probe support costs through comprehensive, cross-functional solutions. The program is built around nearly 40 years of data from proprietary research and over 40,000 successful TEE probe repairs.

TEE Probe Preventive Maintenance

If thousands of devices are covered under your PM plan, shouldn't TEE probes be too?

- Bending rubber failures and deterioration of lens and other seals are the leading cause of costly failures
- Wear to mechanical components leads to reduced performance and articulation failures
- Address points of wear at a **low cost** to minimize risk of costly catastrophic failures

End-User Education

Our clinical experts have crafted educational sessions that provide valuable skills for end-users

- End-users are shown how to identify early warning signs that lead to costly failures
- Service personnel are alerted before gross-fluid invasion occurs

Tools For Success

We've developed processes and tools proven to impact your bottom line

- Interactive forms and webpages guide you through identifying opportunities for improvement
- Custom reports provide actions for clinical areas proven to affect the bottom line and patient throughput

Industry-Leading Capabilities

We're here to assist when needs arise

- Experience and expertise that allows us to confidently say that we have the highest capabilities in the probe repair industry

Repair vs. replacement

Restore safety and performance

Reset your device's lifecycle

Let's rescue your budget, today! Call 844-687-5100 or email safetee@innovatusimaging.com



innovatusimaging.com | 844-687-5100
info@innovatusimaging.com

CENTERS OF EXCELLENCE
PITTSBURGH | TULSA | DENVER

INNOVATUS
imaging
SEE BEYOND